I GENERAL INSTRUCTIONS


1. MANUFACTURER INFORMATION – INTENDED USE

All the products have been manufactured in the facilities of NOVAMIND IKE. Apostolou Pavlou 4 Acharnai Greece.

The website of the company is www.novamind.gr The contact telephone number is +302107770730.

The NOVAMIND dental abutment is intended for use to support single and multiple tooth prosthesis, in the mandible or maxilla. The prosthesis to which the NOVAMIND abutment will connect will be screw retained to the abutment with an anti-rotational internal engagement. The utter aim of the process is for the endosseous implant to be surgically placed in the bone of the upper or lower jaw arches to provide support for the prosthetic devices such as an artificial tooth in order to completely restore patient aesthetics and chewing function.

2. PRODUCT INFORMATION NARRATIVE

3. STORAGE AND HANDLING

All products manufactured by Novamind IKE should be stored at a temperature between 15-25 ° C and 40-60% humidity. The products must be kept away from direct sunlight and any artificial ultraviolet light. The product is well packaged and sealed. A defect on the packaging may involve the loss of the properties of decontamination and disinfection, it is recommended to avoid their use. The material must not be unpacked and handled if it is not going to be immediately used. Implant Novamind IKE products are not sterile. It is recommended sterilization before use, following the methods proposed in section 7.

4. WARNING AND CONTRAINDICATIONS

The reuse of single-use products carries a possible deterioration of its characteristics, which implies the risk of infection of the tissues and /or deterioration of the patient's health. There is a contraindication to the use of the products for patients with conditions that rule out the use of surgery for placement of dental implants. Check the integrity of the packaging, and do not use in case of alteration. NOVAMIND abutments are fully compatible with selected industry dental implants and dental analogues as indicated in NOVAMIND brochures and media. The correct size of the abutment-implant interface is governed by the platform size. Both must have the same platform size so that the abutment and implant fit together. The platform diameter is specified separately on the packaging and is color-coded. If necessary, please check compatibility with your local distributor.

All components of the NOVAMIND abutment series must be protected against aspiration.

The NOVAMIND abutments have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of NOVAMIND abutments in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Use of abutments should be limited to patients with completed craniofacial growth, i.e. patients 18 years old and above without upper age limit exclusions. Deviations to this restriction –if ever- are not endorsed by NOVAMIND, fall exclusively upon the responsibility of the clinician and should be clinically justified on outcome/risk/benefit ratios.

5. RISK ASSOCIATED TO THE USE OF THE PRODUCTS

There is a risk of aspiration or ingestion of the products when used intra-orally, so that appropriate measures be taken to prevent it. If the NOVAMIND abutment or the underlying compatible implant is loaded beyond its functional capacity, excessive bone loss or breakage of the implant or restoration may occur. Clinicians should closely monitor patients for any of the following conditions:

- a. Peri-implant bone loss,
- b. Changes to implant's response to percussion, or
- c. Radiographic changes in bone to implant contact along the implant's length.

Allergic reactions

All materials used in the NOVAMIND Abutment range are biocompatible; however, some patients may present allergies or hypersensitivity to any

GENERAL INSTRUCTIONS FOR USE

Code: NV1a
Approval: 16/05/2018
Version 1
page 2 of 4

of the materials and its components. The abutments and/or auxiliary parts may not be used on patients who are known to have allergies to the corresponding materials. For specific information regarding alloy composition please refer to Section 2.4 herein.

Possible complications

A stressed loading of the implant or abutment over and above its functional capacity can lead to excessive bone loss or fracture of the implant or restoration. The clinician must supervise the occlusion and functional loading of the prosthetic supraconstruction very carefully.

6. CAUTION

Items supplied are intended to be used by qualified health professionals (dental technicians, doctors and dentists). The safety and efficacy of the products supplied by Novamind IKE, is guaranteed only when they are used by trained professionals. Restorations with cantilevers to individual implants are not recommended. Individual restorations with angled abutments should not be used in regions with high mechanical stress. Highly angled abutments on small diameter implants are not recommended for the posterior region of the mouth. Abutments with 3.5 mm diameter are for lateral incisors, mandibular central incisors, or splinting with larger implants. Angled abutments on small diameter implants are not recommended for use in the molar region, unless they are splinted to larger diameter implants. For implants with a small diameter (PF 3.5), the prosthetic restoration should be constructed in such a way that large bending torque does not occur.

7. STERILIZATION METHODS

All NOVAMIND Abutment range products are supplied NON STERILE. For sterilisation, the manufacturer recommends autoclaving the product at 121°C for 30 minutes, drying time 30 minutes (in accordance with standard UNE-EN ISO 17665-1:2007). Most devices are marked for "Single use only" because it is difficult or impossible to clean and decontaminate a used device, reuse can lead to cross-infection. Furthermore, any attempt to reuse a device greatly increases the risk of mechanical failure caused by material fatigue. Any warranty claim resulting from the reuse of a single-use device thus is not accepted by the manufacturer.

Pressure	Temperature [°C]			
	Full Air Free	Removal of 2/3 pf air	Removal of 1/2 of air	No air removal
1/3	109	100	90	72
2/3	115	109	100	90
1	121	115	109	100
4/3	126	121	115	109
5/3	130	126	121	115
2	133	130	126	121

Tabla 1: Influence of incomplete air discharge on autoclave temperature

The autoclave is the device most commonly used for temperatures above 100 ° C. A temperature of 121 ° C (atmosphere pressure) with an exposure time greater than 15 minutes serves to destroy spore forming organisms.

7.1. Advantages

- Fast heating and penetration
- Destruction of bacteria and spores in a short time
- Does not leave toxic waste
- Low deterioration of exposed material
- Economical



7.2. Disadvantages

- It does not allow to sterilize solutions that form emulsions with water
- It is corrosive to certain metallic instruments

GENERAL INSTRUCTIONS FOR USE

Code: NV1a
Approval: 16/05/2018
Version 1
page 3 of 4

8. LABELING SYMBOLS: ISO 15223-1:2012

ISO 15223 Symbology	Description	Novamind Symbology	Description
	Manufacturing date		Recommended torque (N • cm) for screw
	Manufacturer		Torx
	Expiring date		Unigrip
	Do not reuse		Hex (Includes the hex size)
	Batch number		Square drive
	Read instructions for use		Cross drive
	Check instructions for use		TPA Screw
	CE marking, including the number of notified body		
	Product Reference		
	Non Sterile Product		

II SPECIFIC INSTRUCTION FOR USE

1. CASTABLES

Made of Polycarbonate plastic (no waste). Observe the following precautions:

- Use wax to create a sufficient gap around the abutment in order to compensate the higher expansion coefficient of the castable material.
- Screw smoothly to avoid damage.
- Create shapes favoring the filling of cavities and prevent bubbling.
- Use materials for casting with high fluidity.

2. ANALOG

Made of stainless steel AISI-303. For use as a conventional analog must secure a sufficient fixing parts anti-twist and anti-stripping. Verify the compatibility of the connection, in type and size, between the analogue and the prosthetic element prior to tightening. Do not reuse.

3. HEALING ABUTMENT

Made of Titanium grade V. We recommend a gentle hand tightening. The height of this element will be selected to ensures its correct function and prevents transmission of stresses.

4. INTERFACE ABUTMENT

Verify the compatibility of the connection, in type and size, between the interface and the implant. Damage in the connection area of the implant must be avoided. An X-ray in the perpendicular axis to the union interface-implant is recommended in order to ensure the correct adjustment.

5. TITANIUM ABUTMENT

Verify the compatibility of the connection, in type and size, between the titanium abutment and the implant. Damage in the connection area of the implant must be avoided. An X-ray in the perpendicular axis to the union interface-implant is recommended in order to ensure the correct adjustment.

6. Co-Cr CASTABLE ABUTMENT

You need to check compatibility with the implant model you are joining. It should be avoided to damage its area of connection to the implant in the case of carving or machining. The ceramic used with this alloy should have a coefficient of expansion of 14.1×10^{-6} at 500°C or approximate. A pottery far below could give a poor fit, and far above, a break of the same.

7. IMPRESION COPING

Are provided for use in open or closed trays, with or without their specific screws. Prior to use, ensure cleanliness of the implant connection seat. Any dirt could affect the subsequent alignment of the prosthesis. Verify the compatibility of the connection, in type and size, between the impression coping and the implant.

8. SCREWS

Made of Titanium. For fixing abutments on prosthetic implant. It is imperative the strict compliance with the following conditions:

- For screwing or unscrewing use the appropriate screwdriver.
- The screwdriver must be positioned on the longitudinal axis of the joint prosthesis / implant.
- Both for the first fitting a prosthesis, as for future revisions, new screws should be used.
- For immediate load prosthesis:
 - Hand-tight and avoiding excessive torque
 - Prevent rotation of the implant during this operation
- Do not re-use screws from the dental laboratory for clinical use.
- Be sure to use the appropriate reference for each case.
- The minimum number of turns to ensure a good anchorage between five or six, in case of lower number of turns a longer screw shall be used (for more information contact: info@novamind.gr)
- The tightening torque recommended for definitive prosthesis, if not shown on the product label is on the table below:

GENERAL INSTRUCTIONS FOR USE

Code: NV1a
Approval: 16/05/2018
Version 1
page 5 of 4

Model	Torque
01 Group	30 N.cm
02 Group	20 N.cm
03 Group	20 N.cm
04 Group	20 N.cm
05 Group	15 N.cm
06 Group	30N.cm
07 Group	30N.cm
08 Group	30N.cm
10 Group	30N.cm
11 Group	30N.cm
12 Group	30N.cm
33 Group	30N.cm

Table 2, Torque to model recommendation

9. OVERDENTURE ABUTMENT

It is necessary to verify the compatibility with the model of implant to be used. We recommend a radiograph once screwed, at the height of the union with the implant and in the axis perpendicular to that joint, to check the good coupling of the assembly. It is advisable to position the patient so that, if the screw falls during the screwing / unscrewing maneuver, this fall occurs over an area where it can be recovered.

10. RETENTION INSERTS

The retainers must be handled with the corresponding tools distributed by 2018 in order to avoid damage and deformation of the nylon inserts. It is necessary to replace the gums every 6 months or at the time when the patient feels little retention of his prosthesis for overdenture.

NON-STERILIZED PRODUCT: 2018products are not sterilized, but the possibility of sterilizing them is guaranteed

USEFUL INFORMATION ABOUT PRODUCTS MANUFACTURED BY NOVAMIND IKE

**Please read with attention the instructions below to ensure a safe and efficient
use of the products supplied by Novamind IKE**

The entire product range is designed so as to facilitate the work, both in the clinic and dental laboratory, while providing the best quality. Product features and compatibilities are expanded in the commercial catalog; if in doubt please get in touch with us (info@novamind.gr).

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