

Petrol Line.33

Prosthetic components compatible with
Brannemark/3I External Hex Standard platform

33.100

33.



Prosthetic components compatible with
Brannemark/3i External Hex Standard Platform

Lab Analogs.

Ref. Number Diameter \varnothing

33.410.200 4.1




* Only one lab analog for both methods (impression or scanning)

Impression Posts.

Ref. Number Diameter \varnothing Impression Method

33.410.500 4.1 Open Tray



◆ Screw Ref. Number: Open Tray: 33.500.100

Scan Bodies.

Ref. Number Diameter \varnothing

33.410.SA 4.1




◆ Screw Ref. Number: 33.100

* Novamind scan bodies are used only in combination to novamind free libraries either with novamind ti-bases or directly on implant restorative platforms. Scan Bodies are made of peek and are easily scanned with intra-oral and desktop scanners.

Prosthetic components compatible with
Brannemark/3i External Hex Standard Platform

Ti-Bases.

Ref. Number Diameter \varnothing Gingiva Height (mm)



	33.410.600.10	4.1	1.0		Anti-Rotation
	33.410.600.20	4.1	2.0		Anti-Rotation
	33.410.601.10	4.1	1.0		No Anti-rotation
	33.410.601.20	4.1	2.0		No Anti-rotation

Screw Ref. Number: 33.100

* Novamind TI-BASES are colored with a warm Golden Hue for better aesthetic effect. They are available in two types, with and without antirotation depending on the type of restoration (crown or bridge)

Castable Abutments.

Ref. Number Diameter \varnothing



	33.410.300	4.1		Anti-Rotation
	33.410.301	4.1		Anti-Rotation
	33.410.710	4.1		No Anti-rotation
	33.410.711	4.1		No Anti-rotation

Screw Ref. Number: 33.100

* Pre-fabricated Co-Cr implant connection for cast-on.

Titanium Temporary Cylinders

Ref. Number Diameter \varnothing



	33.410.900	4.1		Anti-Rotation
	33.410.901	4.1		No Anti-rotation

Screw Ref. Number: 33.100

* Chairside use or indirect in lab. With and without antirotation depending on the type of restoration (crown or bridge). Golden warm hue for better aesthetics.